Specimen Collection and Submission Guidance for

Molecular Analysis of Confirmed and Suspected Cases of Norovirus

Molecular Tracking of Norovirus to Better Understand Outbreaks and Transmission

Norovirus Stool Specimen Collection

Use: Sterile specimen container with tight-fitting or screw-top lid

Required Specimen: Raw, fresh stool

Required Volume: Minimum of 400 μL, but more is always

preferred

Required Storage and Shipping Temperature:

Store at 4°C. Ship overnight in sealed, insulated containers with cold packs.

- Do not use dry ice.
- Do not freeze specimen.
- Do not use specimen collection kits with preservatives as they interfere with testing.

Norovirus Shipping and Labeling Requirements

Ship as: Category B Biological Substance, UN3373 Specimen must be:

- Triple Packaged to withstand shock, pressure changes, leaks, and other ordinary handling conditions while in transit.
- Packaged with Enough Absorbent Material such as cellulose wadding or paper towels that can soak up the entire contents of the specimen container.
- Ensure all containers are securely closed to prevent leaks!
- Secure lids shut by wrapping in paraffin film (e.g., Parafilm).

Visit DSHS' online <u>Specimen Shipping and Mailing Guidance</u> for more details on shipping Category B substances to the Laboratory.



For specimen size comparison, a single garden pea is approximately 200 μL in volume. Image source: pixabay.com





BIOLOGICAL SUBSTANCE CATEGORY B

Specimen Collection and Submission Guidance for

Molecular Analysis of Confirmed and Suspected Cases of Norovirus

Norovirus Specimens Must be Labeled and Submitted with a G-2B Submission Form

Label Specimen With Unique Identifiers

Every specimen must have at least two unique patient identifiers on its label.



- 2 DOB: 02/19/1993
- 3 06161858


Three patient identifiers provided on this label:

- 1. Name
- 2. Date of Birth
- 3. Medical Record Number

Provide Patient Identifiers in Sections 2 and 3 of Form G-2B

Patient identifiers on specimen label and G-2B submission form must match.

Date of Collection must be provided in Section 3.



	SECTION 3. SPECIMEN				
NOTE: If the 'Date of Collection' field is not completed, the specimen will be rejected.					
	Date of Collection (mm/dd/yyyy) ** 12/21/2023	Time of Collecti	ion ** □ AM √ PM	Col	
G:	Unique Identification Number ** e.g., MRN / Alien # / Accession ID 06161858		Comments or Addition e.g., CDC ID, Previous DSH.		

Select Specimen Type in Section 3

Check "Feces/ stool" as the

specimen type.

Feces / stool Gastric (Aspirate)

Select Test Type in Section 4.4

Check "Norovirus" under *Molecular Studies*.

Select "IDEAS (1610)" in Section 6

Check the IDEAS (1610) box as the Payor.

	☐ Medicaid (2)	☐ Medicare (8)
	Medicaid/Medicare #:	
REQUIRED		
5	☐ Submitter (3)	☐ Immunizations (1609)
g	☐ BIDS (1720)	☐ Private Insurance* (4)
	☐ BT Grant (1719)	☐ TIPP (5144)
*	☐/HIV / STD (1608)	☐ Zoonosis (1620)
	IDEAS (1610)	☐ Other:

Questions About...

Specimen Collection/Suitability: (512) 776-6510

Specimen Shipping: (512) 776-7598 or 1-888-963-7111 ext. 7578 (toll free)

Norovirus Surveillance Program: FoodborneTexas@dshs.texas.gov

Submitter Accounts, Submission Forms, or Result Reports: (512) 776-7578 or LabInfo@dshs.texas.gov

