

INDIVIDUAL SWEEP-CHECK SCREENING (M-40)

Student's name: _____ Age/Grade: _____

School: _____ Teacher: _____

Sweep-Check Screening

1. Instruct and condition each child appropriately for age/grade.
2. Screen three frequencies at 25 dB; begin screening at 1000 Hz.
3. Identify responses with a "+"; identify no response with a "-."
4. Sequence of tone presentations is numbered 1-3 below.

First Screen	Ear	1	2	3	Results
		1000 Hz	2000 Hz	4000 Hz	
	R				☐ Pass ☐ Rescreen ☐ Fail/Refer
Date:	L				

Comments: _____

Signature of Screener: _____ Print Name _____

Children failing to respond to **ONE** (of the three) frequencies in **EITHER EAR** should be rescreened with another sweep-check within three to four weeks. (Signs or symptoms alone would be sufficient for referral.) Failure of **ONE** frequency in either ear on the second sweep-check screen requires a referral or an **Extended Recheck**. If a failure of one frequency occurs when performing the extended recheck, a referral is required.

Second Screen	Ear	1	2	3	Results
		1000 Hz	2000 Hz	4000 Hz	
	R				☐ Pass ☐ Fail
Date:	L				

Comments: _____

Signature of Screener: _____ Print Name _____

Extended Recheck Results

For each of the three frequencies listed, starting at 40 dB; begin screening at 1000 Hz. record the lowest level in decibels (dB) at which the child responds. Record the findings for both the right and left ears. A child should be referred to an appropriately licensed professional if any one of the three frequencies is recorded as greater than 25 dB in either ear.

Extended Recheck	Ear	1	2	3	Results
		1000 Hz	2000 Hz	4000 Hz	
	R	dB	dB	dB	☐ Pass ☐ Fail
Date:	L	dB	dB	dB	

Comments: _____

Signature of Screener: _____ Print Name _____